



Introducing SARS-CoV-2 serology controls

Boost performance at your lab with Thermo Scientific MAS QC solutions

As opposed to kit controls, Thermo Scientific™ MAS™ products are truly independent of any reagents and analyzers used in your lab. This means that the new positive and negative controls for SARS-CoV-2 infections permit **unbiased performance assessment of your serological assays**. With full confidence in your results, you can start optimizing your daily QC and increase your productivity.

#MASproductivity
for COVID-19 testing



21 DAYS
open vial stability @ 2–8 °C



READY TO USE
straight from the fridge



HUMAN
plasma matrix



SARS-CoV-2
Positive Control

Facilitates detection of IgG antibodies
to SARS-CoV-2



SARS-CoV-2
Negative Control

Nonreactive to various infectious agents,
including the SARS-CoV-2 virus




Thermo Scientific MAS Quality Controls for SARS-CoV-2 serology

Part Number	Kit Name	Fill Volume	Quantity per Kit	Shelf Life @ 2-8 °C	Open Vial Stability @ 2-8 °C	Matrix
10028305	MAS SARS-CoV-2 IgG Positive Control Kit	3 mL	5 vials per box	12 months	21 days	Human plasma
10028306	MAS SARS-CoV-2 Negative Control Kit					

Molecular & serological SARS-CoV-2 controls: a powerful combination


Thermo Scientific AcroMetrix™ SARS-CoV-2 controls allow you to verify molecular test performance. Paired with Thermo Scientific MAS controls, you have a complete toolbox for ensuring accurate COVID-19 diagnostics and monitoring in your lab.



Thermo Scientific AcroMetrix SARS-CoV-2 Control

- Mimics human specimen
- For fully automated lab setups

Visit thermofisher.com/sars-cov-2-control



Thermo Scientific AcroMetrix COVID-19 RNA Control

- Full-length genomic RNA
- Low-positive and ultra-low-positive concentrations

Visit thermofisher.com/covid-19-rna-control

Call your local Thermo Fisher Scientific representative or Tech Services for additional information.

Clinical Diagnostics

USA
46500 Kato Road
Fremont, CA 94538 USA
Tel. 800-232-3342 / +1 510-979-5000
sales.diagnostics.fmt@thermofisher.com

Asia Pacific
Tel: +65 6499 9999
cdd.asia.info2@thermofisher.com

Australia
Tel. +61 1800 333 110
auinfo@thermofisher.com

Canada
Tel. +1 905-286-4290
CDD.Canada@thermofisher.com

China
Tel. +86 400 650 5118
sales.china@thermofisher.com

France, Belgium (Brussels / Wallonia) & Luxembourg
Tel. +33 (0) 1 40 86 65 00
cdx.fr.info@thermofisher.com

Germany
Tel. +49 (0) 3302 883 608
cdx.de.info@thermofisher.com

Hong Kong
Tel. +852 3107 7600
SDG.HK.Sales@thermofisher.com

Italy
distributor.cdd@thermofisher.com

Japan
Tel. +81 (0) 120 489 211
JPYOK-CDD.QC@thermofisher.com

Netherlands & Belgium (Flemish)
Tel. +49 3302 883 660
Info.benelux.cdd@thermofisher.com

New Zealand
Tel. +61 0800933 966
nzinfo@thermofisher.co.nz

Nordic
Tel. +46 8 594 691 30
info.nordic.cdd@thermofisher.com

Russia
Tel. +7 495 739 76 41
info.btd.moscow@thermofisher.com

South Africa
Tel. +27 11 792 6790
support-za.idd@thermofisher.com

Spain & Portugal (Iberia)
Tel. +34 93589 8338
CDD.Info.ES@thermofisher.com

Switzerland & Austria
Tel. +41 26 663 86 70
cdx.ch.info@thermofisher.com

United Kingdom & Ireland
Tel. +44 1442 868 940
cdx.UK.info@thermofisher.com

For Countries Not Listed
distributor.cdd@thermofisher.com

Find out more at thermofisher.com/MASproductivity

Availability of product in each country depends on local regulatory marketing authorization status.

© 2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. **CDD-FR-MTL-0249**

ThermoFisher
SCIENTIFIC