



Coming Soon



Pool 5 samples in a single test

Applied Biosystems™ TaqPath™ COVID-19 Pooling Kit



Detect the presence of SARS-CoV-2 in up to five samples simultaneously with the TaqPath COVID-19 Pooling Kit, a real-time PCR multiplex test that enables population testing and pathogen detection using a single kit.

- **Streamline laboratory efficiencies** as the same kit may be used for both pool testing and the deconvolution of positive pools
- **Simplify scaling** your COVID-19 testing by adding to an existing menu and infrastructure
- **Increase processing capacity** by testing up to five samples in a single reaction
- **Employ a complete clinical solution** that offers Applied Biosystems™ COVID-19 Interpretive Software, which converts genetic analysis data into a readable report

Learn more at thermofisher.com/covid19pooling

The content provided herein may relate to products or workflows that have not been officially CE Marked, authorized by the US FDA or global health agencies, and is subject to change without notice; Refer to current user manual for indications, warning and precautions. All validation data to support TaqPath COVID-19 Pooling Kit has been completed and is currently under review by the Food and Drug Administration. Use of TaqPath COVID-19 Pooling Kit is permitted per FDA Policy for Coronavirus Disease-2019 Test during the Public Health Emergency (Revised) released on May 11, 2020. © 2020 Thermo Fisher Scientific Inc. All rights reserved. All trademarks are the property of Thermo Fisher Scientific and its subsidiaries unless otherwise specified. COL013458 1020